

Titanium Micro 540 Needles Derma roller

Usable on the face and body, the 540 Needles Skin Roller may also help cut back and soften pimples or surgical scars, diminish the appearance of cellulite and stretch marks, scale back large pores, clean heavy pitted skin, and assist in reversing the results of sun damage.

The 540 needles Skin Roller dramatically will increase the effectiveness of your hair and pores and skin products by increasing absorption of your products from the usual 4% to as a much as 87% of the product being absorbed deep in your skin. Utilizing the 540 Needles Skin Roller promotes production of missing collagen, resulting in pores and skin rejuvenation, great skin elasticity, and tightening and thickening of the skin. When used on the head (scalp), the Pores and Skin Roller also stimulates the hair to maintain its thickness and help with thinning and balding hair areas. Dramatic results can be seen within one week.

Advantages of 540 Needles Derma roller

- ✓ Rejuvenates skin, promotes pores and skin elasticity, will increase skin thickness.
- ✓ Reduces appearance of cellulite, stretch marks, and large pores.
- ✓ Stimulates hair growth, improves thinning/balding hair areas.
- ✓ Process will not be painful (it's only a tingling sensation).
- ✓ Price effective compared to expensive cosmetic laser treatments.
- ✓ Can be used on all areas of the face, neck and body and scalp.
- ✓ Can be used by yourself anytime at your personal convenience.

How is The 540 Needles Skin Roller Superior to Other Rollers?

Extra Needles per Head The Skin Roller has 540 needles per roller head, in comparison with the standard roller that has less than 200 needles per head. With extra needles, you will have greater penetration in much less time.

Thinner Needle Thickness The needle thickness is 0.20mm, in comparison with the standard needle thickness of 0.25mm, which means you should have much less discomfort throughout needle therapy.

Needle shape leads to less discomfort The diamond shape of the 540 Needles Skin Roller leads to less pain compared to the standard round shaped needle of other skin rollers. This is because the thickness of the spherical-shaped needle is twice as large as the diamond shape needle of the 540 Needles Skin Roller.

Also available in .5 mm and 1.0 mm sizes.

Features and Benefits:

It stimulates dermal healing ability.
{ With no side effects, it's safe and reliable }

The best choice for

- ✓ Anti-aging,
- ✓ Anti Wrinkles,
- ✓ Cell regeneration,
- ✓ Hair regeneration,
- ✓ Acne,
- ✓ Scars,
- ✓ Body Slimming,
- ✓ Cellulite

and much more!

Instructions for use :

Before use, clean your skin, rinse thoroughly and towel dry
Simply use mild pressure gently without forcing the roller into the skin.
Roll 5-6 times per area vertically, horizontally and diagonally
Clean under running warm water after use.
Air-dry.

Important Note :

The roller can be used per one person only!
Keep out of the reach of children

Do not use on

Eczema, Active acne, inflammed, red, oily, acne-prone areas, Psoriasis

Do not share your Roller with anyone else.

Do not use this product if you have open cuts, wounds, sun-burnt skin, active herpes outbreaks, pustule acne lesions, or any other acute infection or inflammation of the skin, have a history of poor wound healing, collagen diseases, blood problems, are pregnant or lactating, are prone to keloid scarring or have diabetes. Aspirin, Nurofen, Vitamin E and blood thinning drugs all cause increased bruising and should not be taken in the two weeks leading up to treatment.

Skin micro-needling (derma rolling) is an increasingly common cosmetic procedure for treating acne scars. Our roller is a class 1 medical device which is simple and easy to use at home. In fact, it has very similar results to Laser Treatment, Chemical Peels and Dermabrasion, but is considerably less invasive (without any of the skin trauma) and carries only a fraction of the cost.

The derma roller perforates the stratum corneum creating a high TEWL (trans epidermal water loss) on acne scars. In a sense it re-activates an inactive healing state allowing ReJuviness scar healing products to work and repair of the stratum corneum. The results will be a totally rebuilt stratum corneum and smooth even colored highly aesthetic re-healed scarred area.

Successful Treatment Regime For Acne scars :

- * Scar Removal Including Pitted and Acne Scarring
- * Boosting the Effectiveness of Scar Healing Products

Suggested Needle Length by Condition:

(.2 - .3 mm) for increased penetration of Hyper Heal scar cream and silicone sheeting into the skin.

(.5 mm) light acne scarring

(1.0 mm) for medium acne scarring

(1.5 mm) for more severe cases of deep acne scars.

(2.0 mm) for professional use in the treatment of acne scarring

Skin Needling Instructions

The instructions are to be used as a guide to getting your skin needling treatment underway. It is important to understand that while a good place to start, it is likely that you will develop your own routine in terms of the frequency and pressure of rolling.

Derma roller Treatments:

The derma roller can be used as part of a successful treatment regime for all of the following.

* Stretch Mark Removal * Anti Aging * Scar Removal including pitted scarring and Acne Scarring * Wrinkle Reduction * Hair Loss Treatments * Boosting topically applied skincare products

There are two recommended methods of using your Derma roller and if you are unsure of how to use them, follow these carefully. We update these recommendations based on the experiences of our customers – whenever the changes occur, we update our web site. Once you have undergone a skin needling session, improvements may continue to develop for several months.

To Start:

Establish a sterile environment, make sure you wash your hands as well as making sure that the body area you are about to treat is also clean. We also recommend that you ensure that the treatment area is clean and free from germs. You may wish to wipe down the area with a disinfectant (rubbing alcohol) or similar.


If it is the first time you are using your Derma roller, sterilizing before use with disinfectant (rubbing alcohol). Although most products do come pre-sterilized you should check before use.

Ensure that you have all the items necessary to begin your treatment including your derma roller, moisturizer, topicalaine gel (if used), antibacterial cream (if used) and any other third party products you may require. (Such as anti cellulite cream, hair loss remedy, etc.)

Method One: The Simple Approach

This is the most common way of using the device - the vast majority of people use their roller in this manner.

Using the roller with slight pressure; roll the device back and forth four times in four directions as shown below. The treatment should not be painful and should not draw blood. You will become accustomed to how much pressure you can apply without causing any pain or discomfort.



Use the roller as part of your daily routine - you can roll every day if you wish although we recommend allowing your skin at least a couple of days of rest each week. It is possible to use the roller far less often and we suggest that fitting your sessions in with your daily and weekly routine will improve the overall experience. Note that it is not necessary to roll more than once a day.

Method Two:

This method involves greater preparation but only requires less frequent use of the device. It is a similar method to that which professional treatments adhere to. Treatments can be carried out whenever you have the opportunity – we recommend monthly treatments although bi-weekly sessions are also fine. When using this method, it is suggested that preparing the treatment area before hand does achieve optimum results. This means applying a scar cream or moisturizer to the area in the

days before treatment will ensure your skin is best prepared to begin the healing process once treatment has finished.

If using a numbing cream such as "EMLA" to the area, please read the instructions on any numbing agents you use thoroughly before use. Follow the instructions above, rolling the Derma roller in different directions, applying greater pressure and/or make additional rolls of the device. (Do not overdo your rolling)

Skin Care After Rolling:

If you have followed method two, we recommend that after cleaning the area, you should apply an antibacterial cream to protect from any possible infection. (It is unlikely this will happen but it is wise to take precautions.) Apply scar cream or moisturizer after use. Remember that part of the benefits of using your roller is the boost in absorption rates of topical products. Nourishing your skin at this point can reap major benefits.

You may experience a slight reddening of the area after treatment. This will fade quickly but dependant on your method, it may last for up to 24 hours. We strongly recommend either staying out of the sun for this period or using a sunscreen to protect yourself.

For Hair Loss and Thinning Hair:

Unlike other treatments, use on thinning hair will not be of benefit without the application of additional products. Use these products as directed, ensuring that you use a short needled derma roller before application (0.5mm or below)

For example - Each morning / night, after showering and washing your scalp, dry well with a clean towel. Roll the Derma roller over the scalp in one direction (and not back and forth). Apply Serum to balding areas and massage into the scalp thoroughly. If you notice the scalp flaking after a few days, simply rub the dead skin away while shampooing your hair. Treat no more than once per day, 5 times per week.

Face Treatments:

Please take care when using the Derma roller on your face and never use the Derma roller on your eyes or lips.

Care of your Derma Roller:

Unless you are disposing your roller after use, when you have finished using the device, rinse it under hot running water. If you have used the roller with numbing cream, stand head first in the half full glass of disinfectant (rubbing alcohol) for approximately one minute. If you only used the roller gently, you need only disinfect the device on a weekly basis. Shake off the excess liquid and stand on a clean towel to dry. Please do not put the device away while it is still wet. When it is completely dry, put the roller away where it will avoid damage.

Important Safety Info:

- * Never share your Dermaroller with other people.
- * Do not use a dermaroller on areas of sunburn, open cuts or other sore and tender areas.
- * Do not use your dermaroller if you are suffering from active herpes outbreaks or any other inflammations of the skin.
- * You can contact us at any time if you have any questions.

* We are not responsible for any unfavourable results which may be derived from improper use of this product.

Possible Side Effects:

Your skin is likely to redden after use. This should fade within a few hours. In extreme cases, this may last for 24 hours.

It is possible that you will experience dryness of skin after rolling. This can be accompanied with some flaking of the skin. If you do suffer in this way, stopping rolling should allow the skin to return to normal within a few days. Please use common sense in these cases - reduce the frequency with which you use the roller and ensure you are using a scar cream or moisturizer after use. If still in doubt or the flaking continues, either stop using the roller or contact us for further support.